

Method 1 – Storage, Retrieval and Transport of biopsies under *Hypothermic* conditions

- a) AQIX® RS-I solution is packaged in 125 mL or 250 mL [PET or glass] containers.
- b) Store AQIX® RS-I solution containers in a fridge @ 3 – 8 °C for not longer than 6 weeks*.
- c) Dispatch AQIX® RS-I solution containers to tissue retrieval site.
- d) Maintain AQIX® RS-I solution containers @ 3 – 8 °C at tissue retrieval site.
- e) Open AQIX® RS-I container for the minimum time possible before inserting the colon biopsy sample (Sample size in each 125 mL container preferably should not exceed 5.0 gm).
- f) Quickly attach the lid closure on the container and close tightly.
- g) Transport the 125 mL container back to laboratory @ 0-4 °C (over ‘wet’ ice) in a polystyrene, outer box.
- h) Conduct experimental procedures immediately upon arrival of the tissue biopsy samples in the laboratory [**recommended**],
or,
- i) Store AQIX® RS-I container + specimen in a fridge @ 3 – 8 °C overnight prior to conducting experimental procedures the next day [**optional**],
or,
- j) Continuously perfuse the isolated colon biopsies in conventional Organ bath equipment using carbogenated AQIX® RS-I solution over a temperature range of 20 – 37 °C

Addendum

1. If tissue contamination at the site of procurement becomes a problem then additional AQIX® RS-I solution may be needed to thoroughly rinse the colon biopsies before inserting into solution specimen bottles. Additionally, 100 mg/L of Chloromycetin may be added to AQIX® RS-I solution to prevent bacterial contamination without compromising the viability of the biopsy specimens.

[* The current working shelf-life for the ‘ready-to-use’ AQIX® RS-I solution is expected to be extended following stability tests now being conducted]